



**YOUR HEALTH
WE CARE !**

MEDTECS GROUP

www.medtecs.com



Medtecs Group is an integrated workwear company and hospital services provider that specializes in the manufacturing and supply of garment products.

Listed on the Singapore Stock Exchange in 1999, the Group's Taiwan Depository Receipts (TDR) was listed on the Taiwan Stock Exchange since 2002.

Group operations commenced in the Philippines in 1989, which have now expanded to China and Cambodia.

Currently Medtecs has 5000 workers across Singapore, Taiwan, Philippines, Cambodia and China.

We have an integrated production line starting from Spinning, Weaving, Dyeing, Printing to Garment for reusable garment. Cleanroom is also available for disposable hospital garment.

Please scan the above QR code for more information.



CAMBODIA

Workforce : 4000
Facilities : Spinning, Weaving,
Dyeing, Printing, Garment



PHILIPPINES

Workforce : 800
Facilities : Weaving, Dyeing,
Garment, Laboratory



CHINA

Facilities : Medical Equipment,
Tray & Kits, Catheters, Bag w/
Tubing, Incontinent Underpad



MANHATTAN SPECIAL ECONOMIC ZON

Service : Land/Facilities
leasing, Labor training &
Management

Reusable Medical Garment 1

Scrub Suit	2
Nurse Uniform	6
Patient Gown	8
Isolation Gown	10
Laundry Bag	10
Infant Shirt	10
Lab Coat	11
Pajama	12
Underpad	14

Personal Protective Equipment 16

Disposable Isolation Gown	17
Disposable Dental Gown	18
Disposable Lab Coat	18
Disposable Warm Up Jacket	19
Disposable Scrub Suit	19
Disposable Bouffant Caps	20
Disposable Face Mask	20
Disposable Coverall	21
Disposable Shoe Cover	21
Disposable Sheet	22

Linen & Blanket 23

Blanket	24
OR Sheet	28

REUSABLE MEDICAL GARMENT





Scrub Suit

V-Neck scrub top
Pants with color-coded drawstring
Reversible for both top and pants

Fiber content: 100% Polyester,
65% Polyester / 35% Cotton,
55% Cotton / 45% Polyester
Fabric weight: 125gsm – 180 gsm
Colors: White and Colors
Size: XS – 5XL



Scrub Suit

V-Neck scrub top
Pants with color-coded drawstring
Reversible for both top and pants

Fiber content: 100% Polyester,
65% Polyester / 35% Cotton,
55% Cotton / 45% Polyester
Fabric weight: 125gsm – 180gsm
Colors: White and Colors
Size: XS – 5XL



Scrub Suit

Fashion style scrub top
Pants with color-coded drawstring

Fiber content: 100% Polyester,
65% Polyester / 35% Cotton,
55% Cotton / 45% Polyester
Fabric weight: 125gsm – 180gsm
Colors: White and Colors
Size: XS – 5XL

Scrub Suit

V-Neck / U-Neck W/buttons/
Fashion style scrub top
Pants with color-coded drawstring
or elastic waist

Fiber content: 100% Polyester,
65% Polyester / 35% Cotton,
55% Cotton / 45% Polyester
Fabric weight: 125gsm – 180 gsm
Colors: White and Colors
Size: XS – 5XL





Nurse Uniform

Fashion styles, high quality, and comfortable.

Fiber content: 100% Polyester,
65% Polyester / 35% Cotton,
55% Cotton / 45% Polyester
Fabric weight: 180gsm – 250gsm
Colors: White and Colors
Size: XS – 5XL





Regular Patient Gown

Complete coverage and modesty
100% Polyester twill tape adds durability
Side-Tie Back, Side-Overlap Back, Straight Back

Fiber content: 100% Polyester,
55% Cotton / 45% Polyester
Fabric weight: 108gsm – 140gsm
Colors: Solid colors or Printed



Regular Patient Gown

Complete coverage and modesty
100% Polyester twill tape adds durability
Butterfly sleeve
Side-Tie Back, Side-Overlap Back, Straight Back

Fiber content: 100% Polyester,
55% Cotton / 45% Polyester
Fabric weight: 108gsm – 140gsm
Colors: Solid colors or Printed



Nursing Gown

Functional gown
2 nursing slits formed by overlapping fabric pieces
4 snaps on I. V. sleeves
Heavyweight construction

Fiber content: 100% Polyester
55% Cotton / 45% Polyester
Fabric weight: 140gsm – 160gsm
Colors: Printed





Isolation Gown

Knit cuff
Tie closure to secure overlap
and prevent gaping

Fiber content: 100% Polyester, 100% Cotton
65% Polyester / 35% Cotton,
55% Cotton / 45% Polyester
Fabric weight: 125gsm – 180gsm
Colors: Yellow and Colors
Size: Free size

Laundry Bag

Buckle closure
Traditional laundry bag style

Fiber content: 100% Polyester
Fabric weight: 183gsm
Colors: Yellow and Colors



Infant Shirt

Fiber content: 100% Cotton,
55% Cotton / 45% Polyester
Fabric weight: 195-220gsm
Colors: White and Colors



Lab Coat

Lapel collar
Snap front closure or with polyester buttons
W/ or w/o Knit cuffs

Fiber content: 100% Polyester, 100% Cotton
65% Polyester / 35% Cotton,
55% Cotton / 45% Polyester
Fabric weight: 180–280gsm
Colors: White and Colors
Size: XS – 5XL



Pajama

Flame-Retardant to BS5722
With lapel collar or V-neck
Long sleeves/Short sleeves

Fiber content: 100% Polyester FR,
55% Cotton / 45% Polyester
Fabric weight: 125gsm
Colors: Solid colors or printed
Size: Large / XL





Nightdress

Flame-Retardant to BS5722

Plain V-neck

Open back /close back style


Fiber content: 100% Polyester FR,

55% Cotton / 45% Polyester

Fabric weight: 125gsm

Colors: Solid colors or printed

Size: Free Size





Underpad

Fiber content: 100% Polyester
Fabric weight: 460gsm, 400gsm
Colors: Blue

Underpad

Surface : 80% Cotton, 20% Polyester, 142gsm

Soaker : 70% Polyester, 30% Rayon, 270gsm

Barrier : 100% PVC, 285gsm

Colors: Orange/Pink/Blue



PERSONAL PROTECTIVE EQUIPMENT





High Risk Isolation Gown with tape

PPSB Coated PE materials provide excellent protection from liquid-borne microbe.
Ultra-sonic welded or heat sealed seams with tape design provides high performance on hydrostatic pressure and fluids impact.
Tie at the neck and waist for secure protection
AAMI level 4 standard

Material: 55% PP / 45% PE
Fabric weight: 36gsm
Colors: Yellow
Size: Free size

Disposable Isolation Gown

Breathable SMS provides fluid-resistant function.
Tie at the waist.
AAMI level 2 standard.

Material: 100% PP SMS
Fabric weight: 25gsm
Colors: Yellow
Style No.: IL-2025-BK



Disposable Surgical Gown

Four-layer SMMS materials, fluid-resistant and breathable
Featured in alcohol-repellency and anti-static treatments.
Hydrostatic Pressure Test reaches 50cm and above.
Ultra-sonic welded and heat sealed seam.
Ties at neck and waist for secure protection
Comfortable knit cuff
Packaged 1ea / bg, 50ea / cs

Fiber content: 100% PP
Fabric weight: 50gsm
Colors: Yellow/Blue
Size: XS – 3XL



Disposable Dental Gown

Made of breathable and fluid-resistant SMS material.
Tie at the waist for wearer to adjust fit.
Designed with knitted cuff.

Material: 100% PP SMS
Fabric weight: 25gsm
Color: White/Yellow
Size: Free size



Disposable Lab Coat

Made of SMS material, breathable,
fluid-repellent and anti-static.
Designed with traditional collar and
knit cuff.

Material: 100% PP SMS
Fabric weight: 35gsm
Colors: Light Blue
Size: M - 4XL



Disposable Warm Up Jacket

Made of soft non-woven material
Knit white cuffs
Front snap closure

Material: 100% PP
Fabric weight: 45gsm-60gsm
Colors: Blue/Green
Size: XS - 5XL



Disposable Scrub Suit

Made of breathable and fluid-resistant
SMS material.
Special treatment for antistatic function.
With one chest pocket and two pockets
at waist.

Material: 100% PP SMS
Fabric weight: 35gsm
Colors: Blue
Size: Regular/ X-Large

Disposable Bouffant Caps

Made of breathable, fluid-repellent PPSB
Puffy shape, typically secured around the head with
elastics.
Packaged 100 ea / bx, 1000 ea / cs

Material: 100% PPSB
Fabric weight: 14gsm
Colors: White/Blue/Green
Style No.:
CA1014W21 / 24 / 27 : White, 21" / 24" / 27"
CA1014B21 / 24 / 27 : Blue, 21" / 24" / 27"
CA1014G21 / 24 / 27 : Green, 21" / 24" / 27"



Face Mask

Premium 3-Ply mask
Bacteria Filtration Efficiency reaches 99.0%
Virus Filtration Efficiency reaches 99.0%
Particulate Filtration Efficiency over 80%
Synthetic Blood Penetration Resistance Test passed
Pleated design

Material: 100% PP
Colors: White/Blue/Green
Style No.:
FM140N-1 : 50ea / bx, 1000 / cs
FM140N-2 : 5ea / bg, 50ea / bx, 1000 / cs





Disposable Coverall

PPSB laminated PE
Breathable and impervious
Synthetic blood penetration resistance test passed
Hood design provides overall protection
With elastic wrist and open ankle

Material: PPSB laminated w/PE film
Fabric weight: 60gsm
Colors: White or colors
Size: XS - XXXL



Disposable Shoe Cover

Non-Skid Shoe Cover/
Heavy Duty Shoe Cover/
Fluid-repellent Impervious Boot Cover

Material: 100% PPSB with PE coating
Fabric weight: 35gsm , 50gsm , 70gsm
Colors: Light Blue
Style No.: SC1035-R : Regular Size
SC1035-XL : X-Large Size

Disposable Transfer Sheet

Disposable transfer sheet saves time and expense of laundering
Intended for single patient use
Helps reduce cross contamination

Material: 100% PPSB
Fabric weight: 90gsm
Colors: Navy Blue
Size: Regular, X-Large



Disposable Stretcher Sheet

Single use
Fluid repellent
Fits standard and large stretchers / gurneys
Elastic ends with deep corner pockets
Reduces laundry/linen costs

Material: 100% PPSB
Fabric weight: 50gsm
Colors: White
Size: Regular, X-Large



Disposable Boot Cover

Made of PPSB coated with PE, non-breathable but fluid-impervious
Ultra-sonic welded seams to avoid blood and fluids penetrating through pinholes
Anti-slip treatment on tread provides additional traction
Adjustable strap ultra-sonic welded to boot

Fiber content: 100% PPSB w/PE coating
Fabric weight: 48gsm
Colors: Blue
Style No.:
SC4048WAB-R White , Regular Size
SC4048WAB-XL White , X-Large Size
SC4048WAB-R Blue , Regular Size
SC4048WAB-R Blue , X-Large Size



BED LINEN & BLANKET





Blanket

With Stripes
1" Hem on Top and Bottom

- Fiber content: 55% Cotton/45% Polyester
- Fabric weight: 260gsm +/- 5%
- Colors: Yellow/Raspberry/Blue/Green
- Size: 180x260 cm



Duo-Tone Blanket ▲

- Fiber content: 100% Polyester
- FR Test Standard: BS EN ISO12952
- Fabric Weight: 369gsm +/- 5%
- Colors: Raspberry/Blue/Green
- Size: 168x244cm




Cellular/Thermal Blanket ▼

2 Sides Hem
Inherent FR yarn

- Fiber content: 100% Polyester FR
- FR Test Standard: BS EN ISO12952
- Fabric weight: 231gsm +/- 5%
- Colors: Raspberry/Green/Blue/Yellow
- Size: 165x210cm

2 Sides Hem
Non-FR

- Fiber content: 55% Cotton, 45% Polyester
- Fabric weight: 300gsm +/- 5%
- Colors: Bleached white
- Size: 168x224cm





Heavy Flannel

- Fiber content: 55% Cotton/45% Polyester
- Fabric weight: 271gsm +/- 5%
- Colors: White/Dyed
- Size: 178x269cm



Ether Blanket

- Fiber content: 55% Cotton/45% Polyester
- Fabric weight: 148gsm +/- 5%
- Colors: Unbleached white
- Size: 178x230cm



Diamond Blanket

- Fiber content: 55% Cotton/45% Polyester
- Fabric weight: 263gsm +/- 5%
- Colors: White/Dyed
- Size: 190x264cm



Blue Stripe Blanket

- Fiber content: 55% Cotton/45% Polyester
- Fabric weight: 212gsm +/- 5%
- Colors: White
- Size: 178x224cm



Sateen Lifter

- Fiber content: 100% Cotton
- Fabric weight: 252gsm +/- 5%
- Colors: Unbleached white
- Size: 137x183, 137x206, 137x230cm



Bath Blanket

1. 1.4 lbs. Blanket
2. Unbleached

- Fiber content: 55% Cotton/45% Polyester
- Fabric weight: 160gsm +/- 5%
- Colors: Unbleached white
- Size: 178x230, 168x230 cm



Snagfree Blanket

- Fiber content: 55% Cotton/45% Polyester
- Fabric weight: 360gsm +/- 5%
- Colors: White
- Size: 188x264cm



Thermal Throw

With 4 Side Hem

- Fiber content: 100% Polyester
- FR Test Standard: BS EN ISO12952
- Fabric weight: 200gsm
- Colors: Light Blue/Green/Lilac/Pink
- Size: 180x210cm

With piping all around

- Fiber content: 100% Polyester
- FR Test Standard: BS EN ISO12952
- Fabric weight: 260gsm +/- 5%
- Colors: Light Blue/Green/Lilac/Pink
- Size: 156x190cm



Thermal Blanket

Heavy weight Thermal Throw

- Fiber content: 100% Polyester
- FR Test Standard: BS EN ISO12952
- Fabric weight: 360gsm
- Colors: Blue
- Size: 168x214cm



Sheets & Pillowcase

T140: 114gsm +/- 5%
T180: 105gsm +/- 5%
T250: 125gsm +/- 5%

- Fiber content: 55% Cotton/45% Polyester
- Colors: White
- Size Sheets: 168x292, 168x274, 168x264 cm
- Size for Pillowcase: 107x81, 107x86 cm



OR Sheet

Fiber content: 100% Cotton
Fabric weight: 150 - 160gsm
Colors: Green



www.medtecs.com